

Radel® dental cases



The sterilization tray

Standard base and lid in injected Radel®

Locking system in stainless steel or integrated


Machined insert in Radel®

Stainless steel basin



Personalization by laser etching, silk screening, etc

Silicone strip laser or water cut



Silicone eyelets
Usable Ø: 2.35, 4.00, 6.50 mm
Choice of colors



Writing desk position possible



AIP MEDICAL

Sterilization trays to your specifications

Radel® dental cases

1563



With rotating turrets

With machined rack
in Propylux®

Internal dimensions:
140 mm x 55 mm



1814



With machined insert

With thermoformed insert

Internal dimensions:
175 mm x 129 mm x 45 mm



2217



With thermoformed insert

With machined insert
and silicone strips

Internal dimensions:
200 mm x 135 mm x 32 mm



2818



With thermoformed
or machined insert

Internal dimensions:
272 mm x 158 mm x 50 mm
Exist also in 34 mm height

Machined racks



Examples of Radel®
and Propylux® racks,
tailor-made machined

