

Radel® bases

The sterilization tray

Presentation of the standard range of Radel® bases. Possibility to combine the bases with our standard lids.

A wide range of dimensions:


- 140 mm x 55 mm x 17 mm
- 175 mm x 129 mm x 20 mm
- 200 mm x 135 mm x 32 mm
- 272 mm x 158 mm x 34 mm
- 270 mm x 173 mm x 40 mm
- 490 mm x 200 mm x 60 mm
- 490 mm x 200 mm x 102 mm




Insert support




Injected or thermoformed Radel®



Sterilization hole \varnothing 4, 6 mm



Customization by laser etching



Foot



Silicone strips and machined parts offer modular solutions



Choice of colors
Additional colors available upon demand



AIP MEDICAL

Sterilization trays to your specifications

D 1563



140 mm x 55 mm x 17 mm*
Sterilization holes Ø 4
Stacking possible

D 2818



272 mm x 158 mm x 34 mm*
Sterilization holes Ø 4
Stacking possible

D 1814



175 mm x 129 mm x 20 mm*
Sterilization holes Ø 4
Stacking possible
Flip-cover

D 3020



270 mm x 173 mm x 40 mm*
Sterilization holes Ø 6
Stacking possible

D 2217



200 mm x 135 mm x 32 mm*
Sterilization holes Ø 6

D 5324



490 mm x 200 mm x 102 mm*
Exist also in:
490 mm x 200 mm x 60 mm*
Sterilization holes Ø 6
Stacking possible

* Internal dimensions